Qualità e design firmato EME  
Quality and design by EME

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## How does it work

**HR-TEK** consists of a generator conveying an alternate current and an application circuit composed of two electrodes.

These electrodes are applied on the patient's body that becomes a part of the electric circuit.

The electrodes can be:

- **resistive type electrodes**: the transfer is made through a resistive (non isolated) electrode mobilizing the charges, thus making it possible to concentrate them on deeper and more resistive areas (bone matrix and deep muscles)

- **capacitive type electrodes**: the capacitive contact transfer is made through an electrode covered in an insulating ceramised material, available in different dimensions depending on the area to be treated, mobilizing the charges (ions) in the subcutaneous tissues.

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## The main fields of application are the following:

- Dermatologic pathologies
- Traumatology
- Sports medicine
- Aesthetic medicine
- Skeletal muscle pain
- Inflammatory pathologies of joints and tendons
- Peri-articular pathologies (bursitis, calcifications, etc.)

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**HR-TEK** devices are certified  
KIWA CERMET ITALIA Spa CE 0476, according to Directive 93/42/EEC

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**The tecartherapy**

*The capacitive/resistive endothermy or Tecartherapy, (acronym for capacitive and resistive energetic transfer) is a therapy stimulating the natural reparation processes of the organism and shortening motor recovery time. Diathermy, through the electromagnetic energy-tissue interaction, produces an increase in temperature occurring inside the tissues in an uniform and controlled way. Such electromagnetic interaction causes a ionic flow with a micro-hyperaemia, which, finally, encourages the release of endogen "substances" (most of all cortisol and endorphins) allowing the reduction of pain, edemas and inflammations. Therefore, HR-TEK stimulates an increase in the hematic flow in a direct way thanks to the increase in temperature and in an indirect way through an oxygen demand from the treated tissues. The increase in the hematic flow encourages the augmentation in the normal immune defences and stimulates tissue regeneration. The HR-TEK diathermy is a system optimising such transmission to body tissues in a safe way for the operator and the subject of the treatment in terms of:*


- EFFECTIVENESS
- RAPIDITY OF RESULTS

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**Operating mode and benefits**

HR-TEK operates in two modes:

- 1) Capacitive mode:** the action is concentrated on the area underneath the mobile electrode. It actively acts on soft tissues rich in water content, such as muscles, the venous-lymphatic system, ecc. HR-TEK, when used in capacitive mode, gives excellent results for:
- The treatment of fracture aftereffects
  - Articular and peri-articular pathologies
  - The drain of post-traumatic edemas

- 2) Resistive mode:** it performs an effective action on all higher resistance tissues, poor in water content. The biological effect occurs in the higher resistance tissues interposed between the mobile electrode and the return plate, such as bones, joints, tendons, ligaments, cartilages. HR-TEK, when used in resistive mode, gives excellent results for:
- Generic traumatology
  - Sports traumatology
  - Arthritis
  - Arthrosis
  - Tendinitis
  - Capsulitis and bursitis
  - Calcaneal spur syndrome